customerservice@carolinaflow.com •

www.carolinaflow.com

919.635.8438



FLOWTAINER™ CLOSURE SOLUTIONS

When a single drop of fluid can make the difference in someone's health, you need smart fluid handling solutions that get the job done safely. Manufacturers worldwide rely on connector products that improve yield, cut costs and reduce time to market.



FlowTainer™ closure solutions are ideal for critical fluid handling in upstream and downstream processing, such as transfer, sampling, freezing, thawing, storage, and transport. Closure assemblies can be uniquely customized to solve your process challenge. FlowTainer™ closure solutions are designed to maintain an aseptic seal and closed sterile system, eliminating the need for cleaning validation and ensuring product integrity.

TABLE OF CONTENTS

- GL45 Caps and Assemblies for Glass Containers | Page 2
- 83B Caps and Assemblies for Carboys (Nalgene® and Foxx LS) | Page 4
- 50mL and 15mL Caps and Assemblies for Centrifuge Tubes (Labcon®) | Page 6
- ◆ Ultra Clean Bottles | Page 8

GL45 Caps and Assemblies for Glass Containers

SPECIFICATIONS:

| ~ | Manufactured in an ISO 8 Cleanroom |
|----------|------------------------------------|
| | |

✓ ISO 9001:2015 Certified Facility

USP Class VI and animal component free

Sterilizable by autoclave and gamma irradiation

Assemblies come pre-packaged ready for autoclaving




| PART NUMBER | DESCRIPTION | CAP MATERIAL | HOSE BARB ID | CLOSURE SIZE | PACK SIZE |
|-------------|------------------|--------------|--------------|--------------|-----------|
| F-BCC-45-1 | 2 port Cap, GL45 | PP | (2) - 1/8" | GL45 | 10 |
| F-BCC-45-2 | 2 port Cap, GL45 | PP | (2) - 1/4" | GL45 | 10 |
| F-BCC-45-3 | 3 port Cap, GL45 | PP | (3) - 1/8" | GL45 | 10 |


| PART NUMBER | DESCRIPTION | TUBING | CONFIGURATION | | STERILITY CLAIM | FIGURE DETAILS |
|-----------------|---------------------------|----------------------------|---|------|--------------------|-------------------|
| F-BCCA-45-2-PCS | 2 port Cap Assembly, GL45 | Platinum Cured Silicone | | Each | Non-Sterile | Figure 1 |
| F-BCCA-45-2-TPE | 2 port Cap Assembly, GL45 | AdvantaFlex® TPE | Line 1: 1/8" x 1/4"; 24" OAL w/ 12" dip tube Line 2: 1/8" x 1/4", 4" OAL | Each | Non-Sterile | Figure 1 |
| F-BCCA-45-2-PM | 2 port Cap Assembly, GL45 | PharMed® | | Each | Non-Sterile | Figure 1 |
| F-BCCA-45-3-PCS | 3 port Cap Assembly, GL45 | Platinum Cured Silicone | Line 1: 1/8" x 1/4"; 24" OAL w/ 12" dip tube | Each | Non-Sterile | Figure 2 |
| F-BCCA-45-3-TPE | 3 port Cap Assembly, GL45 | AdvantaFlex® TPE | Line 2: 1/8" x 1/4"; 24" OAL w/ 12" dip tube Line 3: 1/8" x 1/4", 4" OAL | Each | Non-Sterile | Figure 2 |
| F-BCCA-45-3-PM | 3 port Cap Assembly, GL45 | PharMed® | | Each | Non-Sterile | Figure 2 |

Please see Figure Details on the following page.



GL45 Caps and Assemblies for Glass Containers







83B Caps and Assemblies for Carboys (Nalgene® and Foxx LS)

SPECIFICATIONS:

- Manufactured in an ISO 8 Cleanroom
- ISO 9001:2015 Certified Facility
- USP Class VI and animal component free
- Sterilizable by autoclave and gamma irradiation
- Available with sterility assurance level 10⁻⁶ per ISO 11137




| PART NUMBER | DESCRIPTION | CAP MATERIAL | HOSE BARB ID | CLOSURE SIZE | PACK SIZE |
|-------------|-----------------|--------------|------------------------|--------------|-----------|
| F-BCC-83-1 | 3 port Cap, 83B | PP | (2) - 1/2", (1) - 1/4" | 83B | 10 |
| F-BCC-83-2 | 3 port Cap, 83B | PP | (2) - 3/8", (1) - 1/4" | 83B | 10 |
| F-BCC-83-3 | 3 port Cap, 83B | PP | (3) - 1/4" | 83B | 10 |

| PART NUMBER | DESCRIPTION | TUBING | CONFIGURATION | QUANTITY | STERILITY CLAIM | FIGURE DETAILS |
|-----------------|---|----------------------------|---|----------|--------------------|-------------------|
| F-BCCA-83-2-PCS | 3 port Cap, 83B, for 5L Carboy, PCS | Platinum Cured Silicone | Line 1: 3/8" x 5/8"; 24" OAL w/AseptiQuik® G, 12" dip tube | Each | Sterile | Figure 1 |
| F-BCCA-83-2-TPE | 3 port Cap, 83B, for 5L Carboy, TPE | AdvantaFlex® TPE | Line 2: 3/8" x 5/8"; 24" OAL w/AseptiQuik® G, 12" dip tube | Each | Sterile | Figure 1 |
| F-BCCA-83-2-PM | 3 port Cap, 83B, for 5L Carboy, PharMed® | PharMed® | Line 3: 1/4 x 1/2", 4" OAL w/ .2µm vent filter | Each | Sterile | Figure 1 |
| F-BCCA-83-2-PCS | 3 port Cap, 83B, for 10L Carboy, PCS | Platinum Cured Silicone | Line 1: 3/8" x 5/8"; 24" OAL | Each | Sterile | Figure 1 |
| F-BCCA-83-2-TPE | 3 port Cap, 83B, for 10L Carboy, TPE | AdvantaFlex® TPE | w/AseptiQuik® G, 12" dip tube Line 2: 3/8" x 5/8"; 24" OAL w/AseptiQuik® G, 12" dip tube Line 3: 1/4 x 1/2", 4" OAL w/ .2µm vent filter | Each | Sterile | Figure 1 |
| F-BCCA-83-2-PM | 3 port Cap, 83B, for 10L Carboy, PharMed® | PharMed [®] | | Each | Sterile | Figure 1 |
| F-BCCA-83-2-PCS | 3 port Cap, 83B, for 20L Carboy, PCS | Platinum Cured Silicone | Line 1: 3/8" x 5/8"; 24" OAL | Each | Sterile | Figure 1 |
| F-BCCA-83-2-TPE | 3 port Cap, 83B, for 20L Carboy, TPE | AdvantaFlex® TPE | w/AseptiQuik® G, 12" dip tube Line 2: 3/8" x 5/8"; 24" OAL w/AseptiQuik® G, 12" dip tube | Each | Sterile | Figure 1 |
| F-BCCA-83-2-PM | 3 port Cap, 83B, for 20L Carboy, PharMed® | PharMed° | Line 3: 1/4 x 1/2", 4" OAL w/ .2µm vent filter | Each | Sterile | Figure 1 |

Please see Figure Details on the following page.



83B Caps and Assemblies for Carboys (Nalgene® and Foxx LS)





50mL and 15mL Caps and Assemblies for Centrifuge Tubes (Labcon®)

SPECIFICATIONS:

| V | Manufactured in an ISO 8 Cleanroom |
|------------|---|
| ⊘ | ISO 9001:2015 Certified Facility |
| \bigcirc | USP Class VI and animal component free |
| *** | Sterilizable by gamma irradiation |
| 攻 | Aseptically connect using TPE tube welder |
| Ø | Available with sterility assurance level 10-6 per ISO 11137 |


| PART NUMBER | DESCRIPTION | MATERIAL | HOSE BARB ID | CLOSURE SIZE | PACK SIZE |
|-------------|-----------------------------------|----------|--------------|--------------|-----------|
| F-CC-50ML | 50mL Tube Cap, 2 port | HDPE | (2) - 1/8" | 50mL | 50 |
| F-CC-15ML | 15mL Tube Cap, 2 port | HDPE | (2) - 1/8" | 15mL | 50 |
| 00113096 | 15 mL PerformR° Centrifuge Tube | PS | N/A | 15mL | 125 |
| 00113097 | 50 mL SuperClear® Centrifuge Tube | PP | N/A | 50mL | 50 |
| 00114530 | 15 mL SuperClear® Centrifuge Tube | PP | N/A | 15mL | 50 |


| PART NUMBER | SAMPLE VOLUME | CONTAINER MATERIAL | CONFIGURATION | QUANTITY/ PACK | STERILITY Claim | FIGURE Details |
|----------------|------------------|-----------------------|--|-------------------|--------------------|-------------------|
| F-CTA-PP-15D | 15mL | Polypropylene | | 5 | Non-Sterile | Figure 1 |
| F-CTA-PS-15D | 15mL | Polystyrene | Line 1: 1/8" x 1/4"; 12" OAL w/ dip tube | 5 | Non-Sterile | Figure 1 |
| F-CTA-PP-15D-S | 15mL | Polypropylene | Line 2: 1/8" x 1/4", 3" OAL w/ vent filter | 5 | Sterile | Figure 1 |
| F-CTA-PS-15D-S | 15mL | Polystyrene | | 5 | Sterile | Figure 1 |
| F-CTA-PP-15 | 15mL | Polypropylene | | 5 | Non-Sterile | Figure 2 |
| F-CTA-PS-15 | 15mL | Polystyrene | Line 1: 1/8" x 1/4"; 12" OAL | 5 | Non-Sterile | Figure 2 |
| F-CTA-PP-15-S | 15mL | Polypropylene | Line 2: 1/8" x 1/4", 3" OAL w/ vent filter | 5 | Sterile | Figure 2 |
| F-CTA-PS-15-S | 15mL | Polystyrene | | 5 | Sterile | Figure 2 |
| F-CTA-PP-50D | 50mL | Polypropylene | Line 1: 1/8" x 1/4"; 12" OAL w/ dip tube | 5 | Non-Sterile | Figure 3 |
| F-CTA-PP-50D-S | 50mL | Polypropylene | Line 2: 1/8" x 1/4", 3" OAL w/ vent filter | 5 | Sterile | Figure 3 |
| F-CTA-PP-50 | 50mL | Polypropylene | Line 1: 1/8" x 1/4"; 12" OAL | 5 | Non-Sterile | Figure 4 |
| F-CTA-PP-50-S | 50mL | Polypropylene | Line 2: 1/8" x 1/4", 3" OAL w/ vent filter | 5 | Sterile | Figure 4 |


Please see Figure Details on the following page.




50mL and 15mL Caps and Assemblies for Centrifuge Tubes (Labcon®)









ULTRA CLEAN BOTTLES

DESCRIPTION:

FlowTainer™ ultra clean bottles are the optimal solution for critical fluid handling in upstream and downstream processing, such as filling, freezing, thawing, storage, and transport. The robust impact resistant bottles are made from USP Class VI and animal ingredient free polycarbonate. Manufactured in an ISO 7 cleanroom, ultra clean bottles undergo stringent testing to ensure the highest quality standards are met, such as USP <85> and USP <788>.

FlowTainer™ ultra clean bottles are available from 5 mL to 20 L with a variety of hose barb cap configurations for use with platinum cured silicone and TPE tubing. Choose from standard bottles and assemblies or customize to fit your unique process application.

SPECIFICATIONS:



ISO 9001:2015 Certified Facility

Easy to read graduations

Ultra-low temperature storage and transport to -90°C

USP Class VI and animal component free

Sterilizable by gamma irradiation and autoclave

Excellent for freezing and storage applications

Validation package with extractables available

Available with sterility assurance level 10⁻⁶ per ISO 11137





ULTRA-LOW TEMPERATURE STORAGE



ULTRA CLEAN BOTTLES

FLOWTAINER™ REGULATORY OVERVIEW:

Material of Construction Polycarbonate (PC) Bottle and HDPE Cap **Biocompatibility** USP <85>, Bacterial Endotoxin USP <87>, Biological Reactivity Tests, Cytotoxicity, In Vitro USP <88>, Biological Reactivity Tests, Class VI, In Vivo ISO 10993-4, Hemolysis **Physiochemical** USP <661>, Plastic Packaging Systems FDA 21 CFR 177.1580, Polycarbonate Resins FDA 21 CFR 177.1520, Olefin Polymers **Particulate** USP <788>, Particulate Matter in Injections **Extractables** Per guidelines of BPOG and USP <665> Legislation REACH RoHS **Shelf Life**



3 Years, Non-Sterile/Non-Irradiated

ULTRA CLEAN BOTTLES

PART NUMBERS:

| PART NUMBER | VOLUME | BOTTLE MATERIAL | CAP MATERIAL | CLOSURE SIZE | STERILITY CLAIM | QTY PER PACK/CASE |
|---------------|--------|-----------------|--------------|--------------|-----------------|-------------------|
| F-PCB-005ML | 5 mL | PC | HDPE | 24 mm | Non-Sterile | 25/500 |
| F-PCB-005ML-G | 5 mL | PC | HDPE | 24 mm | Sterile | 25/500 |
| F-PCB-020ML | 20 mL | PC | HDPE | 24 mm | Non-Sterile | 25/500 |
| F-PCB-020ML-G | 20 mL | PC | HDPE | 24 mm | Sterile | 25/500 |
| F-PCB-125ML | 125 mL | PC | HDPE | 48 mtm | Non-Sterile | 5/50 |
| F-PCB-125ML-G | 125 mL | PC | HDPE | 48 mm | Sterile | 5/50 |
| F-PCB-500ML | 500 mL | PC | HDPE | 48 mmt | Non-Sterile | 2/48 |
| F-PCB-500ML-G | 500 mL | PC | HDPE | 48 mm | Sterile | 2/48 |
| F-PCB-001L | 1L | PC | HDPE | 48 mm | Non-Sterile | 2/28 |
| F-PCB-001L-G | 1L | PC | HDPE | 48 mm | Sterile | 2/28 |
| F-PCB-002L | 2 L | PC | HDPE | 48 mm | Non-Sterile | 1/20 |
| F-PCB-002L-G | 2 L | PC | HDPE | 48 mm | Sterile | 1/20 |
| F-PCB-005L | 5 L | PC | HDPE | 70 mm | Non-Sterile | 1/6 |
| F-PCB-005L-G | 5 L | PC | HDPE | 70 mm | Sterile | 1/6 |
| F-PCB-010L | 10 L | PC | HDPE | 70 mm | Non-Sterile | 1/3 |
| F-PCB-010L-G | 10 L | PC | HDPE | 70 mm | Sterile | 1/3 |
| F-PCB-020L | 20 L | PC | HDPE | 70 mm | Non-Sterile | 1/3 |
| F-PCB-020L-G | 20 L | PC | HDPE | 70 mm | Sterile | 1/3 |

CLOSURES FOR FLOWTAINER™ ULTRA CLEAN BOTTLES

| PART NUMBER | DESCRIPTION | CAP MATERIAL | HOSE BARB ID | CLOSURE SIZE | STERILTY CLAIM |
|-------------|-------------------|--------------|------------------------|--------------|----------------|
| F-BCC-48-1 | 2 Port Cap, 48 mm | HDPE | (1) - 1/8", (1) - 1/4" | 48 mm | Non-Sterile |
| F-BCC-48-2 | 2 Port Cap, 48 mm | HDPE | (2) - 1/4" | 48 mm | Non-Sterile |
| F-BCC-48-3 | 2 Port Cap, 48 mm | HDPE | (2) - 1/8" | 48 mm | Non-Sterile |
| F-BCC-48-4 | 3 Port Cap, 48 mm | HDPE | (3) - 1/4" | 48 mm | Non-Sterile |
| F-BCC-48-5 | 3 Port Cap, 48 mm | HDPE | (2) - 1/4", (1) - 1/8" | 48 mm | Non-Sterile |
| F-BCC-70-1 | 2 Port Cap, 70 mm | HDPE | (2) - 1/4" | 70 mm | Non-Sterile |
| F-BCC-70-2 | 3 Port Cap, 70 mm | HDPE | (3) - 1/4" | 70 mm | Non-Sterile |
| F-BCC-70-3 | 4 Port Cap, 70 mm | HDPE | (4) - 3/8" | 70 mm | Non-Sterile |
| F-BCC-70-4 | 4 Port Cap, 70 mm | HDPE | (4) - 1/4" | 70 mm | Non-Sterile |

